



**Bay Alarm
Medical**



SOS

In-Home Cellular

Medical Alert System

USER GUIDE



Table of Contents

- 4** Setting Up a Cellular System
- 4** Testing your SOS Home
- 5** Automatic Fall Detection
- 6** Customer Referral Benefit

Emergency dispatchers
may call. Save this
number for caller ID
recognition:

1-844-996-0208

Customer Service & Tech Support

1-877-522-9633

Instructional Videos

www.BayAlarmMedical.com/videos

What's Included

SOS Home Base Kit



In-Home Base Unit
w/AC Power Adapter
Cellular



Personal Help Button
with Wrist band & Neck

What you need to set up



AC Power Outlet
not controlled by a switch

Add on accessories



**Auto Fall
Detection (AFD)**



Wall Buttons



Bella Charms
for standard buttons





Lock Box

Setting Up A Cellular System

Follow all the set-up directions to ensure proper system functions.

1. Place the **Base Unit** in an area with strong signal (i.e. by a window, or a glass door).
2. Plug the Base Unit's power cord into the power outlet (**that is not controlled by a switch**).
3. Find the '**O/I**' power switch on the rear of the Base Unit and slide it to the '**I**' position.
4. The status light will start to flash **orange** while it searches for cellular coverage
5. When ready, the Base Unit will announce, "**System Ready**" and the status light will turn **green**.

Wait 30 seconds before moving on to the next steps.



Testing your SOS Home System

The **Range Test** will help determine the connection between your buttons and the Base Unit in and around your home. Please perform the following **for each of your help buttons** before using.

NOTE: It may be helpful to have two people perform the range test. **Person 1** will stay at the Base Unit to listen for the steady tone. **Person 2** will press the personal Help Button in different locations around the residence and yard.

1. Press the '**T/L**' button, on the rear of the **Base Unit**. It will announce, "**Range Test Mode.**"
Device will time out and reset if no action is taken after 45 seconds.
2. Press the **Personal Help Button**. A **steady tone will sound from the Base Unit**.
This means the unit is detecting the button.

Repeat these steps for various locations throughout the residence and yard.

3. Once areas in and around the residence are tested, press the '**RESET**' button. Unit should announce, "**System Ready.**"

NOTE: If you press the Personal Help Button and the Base Unit **does not follow with a tone**, this means your button is not detectable at that location and **will not trigger a help signal**.

Section Continues to page 6 >

Test the Device

The **Device Test** will help you confirm that your unit is fully functional.

NOTE: We recommend that you perform a Device Test **once a month**.

1. Press the 'HELP' button on the Base Unit. It will announce, **"Calling for help."**
2. After a period of beeps, followed by brief silence, **a live operator will connect over the 2-way speaker** on the Base Unit. Let the operator know that you are testing. Follow their instructions.
3. Once the test is complete, press the **'RESET' button** on the Base Unit.

NOTE: Follow the same steps to test your personal Help Button and any wall button.

NOTE: If you do not hear an operator from the 2-way speaker, you should receive a call on the phone number provided to us.

If you get neither response from the 2-way speaker nor from your home telephone within 3 minutes, call us immediately.



Auto Fall
Detection
Only



Understanding Automatic Fall Detection

Automatic Fall Detection (AFD) is designed to call for help when you're unable to. Loved ones with a history of falling may find this feature beneficial.

This service may not be appropriate for everyone.

How It Works

Auto Fall Detection (AFD) Activation

The device is worn on the outside of one's clothing. AFD is triggered when there's sudden, downward movement, followed immediately by lack of movement.

Lack of movement

Stillness communicates to the sensor that a fall has occurred, that the user is immobilized and can't physically push a button.

If you can, press button

For the best efficiency, always press the button if you're physically able to do so, even if an auto fall alarm is triggered.

Ways to prevent false alarms

Depending on lifestyle, the Automatic Fall Detection (AFD) could be set off unintentionally by movements that mimic a fall. The following will help prevent that from occurring.

Avoid certain quick movements

Avoid fast movements that can trigger a false alarm, such as plopping down on a couch, or leaning over to pick up something.

Do not wear it while sleeping

The AFD button should not be worn to bed. Sleeping motions can set off a false alarm. When this happens, if you don't reply to **the Monitoring Center**, then an emergency response will be sent to the home.

Handle the device with care

Be careful while putting on or taking off the AFD button. If not wearing it, place it in a location where it won't drop or dangle.

If a False Alarm Occurs

Sometimes false alarms can happen even while taking the necessary precautions. If this happens, please let the call go through to **the Monitoring Center**. A representative will respond from **the Base Unit**. Please respond back to let them know that you're okay, that it was a false alarm.



Refer A Friend, Reward Yourself

Receive a **FREE** month of service if you have a friend sign up with Bay Alarm Medical.

FAQs

How do I get my free month?

After 60 days following your friend's activation, we will process a credit for a **free month of service to you**.

In order to receive your credit, please make sure they mention your full name at the time of order.

Is there a referral limit?

Not at all! Tell as many friends as you'd like. Receive a free month **for each friend that signs up with a new subscription plan**.

If 12 friends sign up and stay with us for 60+ days, then you could earn 12 months of free service!



**Bay Alarm
Medical**

SOS

In-Home Cellular

Medical Alert System

www.BayAlarmMedical.com

© All Rights Reserved. Bay Alarm Medical & GetSafe.

202402