



## SOS Smartwatch

Enhance your life with the modern day alternative to traditional medical alert systems. The Bay Alarm Medical SOS Smartwatch's crystal clear touch screen and simple functionality will allow you to make the most of your day while staying protected. Whether your in or outside the home, the SOS Smartwatch's portability will make it easy to wear wherever you go and do.



### Features

- Two-Way Speaker
- Physical SOS Button
- Powered by 4G LTE
- Step Counter
- 3 Screen Options
- No Smartphone Required!
- Weather Proof

### Specifications

Screen Size	Wristband Length	Wristband Width	Weight	Battery of Life
1.25in	5.25-7.25in (8in w/ extension)	0.875in	~2oz	Up to 18-24 hrs

### Specifications

Placing an emergency call with the SOS Smartwatch's physical button is simple to do. During an emergency, press and hold the SOS button on the left hand side for 5 seconds. The watch will connect to our 24/7 emergency center. Once connected, the dispatcher will ask if help is needed and with the two-way speaker, you can communicate your situation. If help is needed, the built-in GPS will pinpoint your location and send necessary help your way.



### Fits Most Wrist Sizes

The SOS Smartwatch is available in two wristband sizes. Before purchasing, please measure your wrist to determine which size will best fit you or your loved ones.

**Size One** = 5.25 - 7.25 inches

**Size Two** = 5.25 - 8.00 inches

Check out our videos or instructions for further information at [bavalarmmedical.com/videos/](http://bavalarmmedical.com/videos/) and [bavalarmmedical.com/help-guides/](http://bavalarmmedical.com/help-guides/) or call 1-877-522-9633 to order.